







ORTOPEDIJOS TECHNIKA

MATERIALS

2025-1



Leathers	1
Textiles	8
Linings	10
Paddings	12
Stiffenings	13
Fastenings	15
Welts	16
Midsole materials	17
Boardings	18
Soles	19
Molded soles	25
Vacuum insoles	26
Milled insoles	28
Insoles top covers	29

Please note:

- colour shades of materials may vary.
- the colors of materials may vary between real life and photographs due to photo quality, screen resolution, and/or print quality if this palette is printed.

LILY

Thickness
1.0-1.2 mm



SETA

Thickness
1.3-1.5 mm



TULIP

Thickness
1.1-1.3 mm




ART

Thickness
1.1-1.2 mm



SNAKY

Thickness
0.9-1.1 mm



MOUSSE

Thickness
1.0-1.1 mm



VELVET

Thickness
0.9-1.0 mm



VOLCANO

Thickness
0.9-1.1 mm



LUMIÈRE - LIPSTICK

Thickness
0.9-1.1 mm




CRONOS

Thickness
0.9-1.1 mm



VANESSA

Thickness
0.8-1.0 mm



PURESOFT

Thickness
1.1-1.3 mm




COCCO OYSTER

Thickness
0.9-1.0 mm



SOUPLESSE

Thickness
0.9-1.0 mm



BOMBAY

Thickness
0.9-1.0 mm



JAIPUR PEARL

Thickness
0.9-1.1 mm



VEDETTE

Thickness
1.1-1.3 mm



KIDEPO

Thickness
0.8-0.9 mm



MEDICI

Thickness
1.2-1.4 mm



BOXLIGHT

Thickness
1.0-1.2 mm



COLORADO

Thickness
1.1-1.3 mm



IMPRUNETÀ

Thickness
1.4-1.6 mm



SPORTLUX

Thickness
1.3-1.5 mm



HUNTER

Thickness
1.6-1.8 mm



SURVIVAL HYDRO

Thickness
1.4-1.6 mm



ESPRIT

Thickness
1.2-1.4 mm



TIFFANY HYDRO

Thickness
1.3-1.5 mm



ONTARIO HYDRO

Thickness
1.5-1.7 mm



PAMPA HYDRO

Thickness
1.7-1.9 mm



HEAVY DUTY

Thickness
2.0-2.2 mm



DIACALF

Thickness
0.8-1.0 mm



CRACKED

Thickness
0.9-1.1 mm



STING

Thickness
1.1-1.3 mm



COINS

Thickness
0.9-1.1 mm



SPECIAL MATERIALS




Waterproof

BREATHABLE


TX21
White




TX20
Light grey




TX22
Gris grey




TX10
Red



TX23
Azul blue




TX24
Marine



TX07
Black




TX42
Grey




ELASTIC


TX27
Chaki with padding




TX25
Jeans with padding




TX26
Black with padding



TX18
Elastic black



TX19
Tricot negro



TX29
Stretchable black





DURABLE

TX30
Camouflage



TX31
Cordura
black



TX32
Cordura
brown




TX33
Cordura
blue




TX34
Cordura
chaki



TX35
Cordura
grey



TX37
Cordura
granate




TX38
Cordura
Teja



RECYCLED

TX28
Recycled
cotton
black



Recycled material

NATURAL LEATHER

LN02

Calflining cream,
0.7-0.9 mm



LN45

Natura lining burro,
1.0-1.1 mm



LN26

Lamnappa naturel,
0.7-0.9 mm



LN22

Calflining camel,
0.9-1.1 mm



LN24

Calflining nero,
0.7-0.9 mm



SOFT

LN43

Foam laminated suede beige,
4 mm



LN38

Foam laminated suede chocolate,
4 mm



LN39

Foam laminated suede nero,
4 mm



THERMO

LN28

Natural fur off white



LN30

Natural fur brown



LN52

White thermo with membrane



LN27

Natural fur black



LN42

Puratex,
3.2 mm



LN44

Foam laminated velour dar grey,
4 mm



Waterproof

BREATHABLE

LN47
DuraFit
PLUS,
1.1 mm



LN06
Foam
laminated
PE flisse
black,
4 mm



TX39
Textile black,
3 mm



RECYCLED

LN46
Recycled
camel,
1 mm



Recycled material



PU FOAM

SFM 04-81
3 mm

SFM 05-81
5 mm

SFM 06-81
8 mm




ZEG

SFM 07-81
3 mm


SFM 08-81
5 mm




SOFT	Erkoflex			STM 20 (1,8 mm)	Erkoflex			STM 10 (3 mm)	TEPP2			STM 21 (1 mm)
	Eva 35 Sh			STM 27 (3 mm)	TALYN			STM 28 (1 mm)				
	Interflex			STM 09 (1 mm)	MAX			STM 15 (1,5 mm)	MAX			STM 16 (1,1 mm)
	Leather			STM 11 (1,8-2 mm)	Rhenoflex			STM 17 (1,5 mm)	Rubberflex			STM 26 (1,7 mm)
	Ibiflex			STM 12 (2 mm)	Thermoflex			STM 14 (1,5 mm)				







HEEL COUNTER




TOE CAP




BRANDSOLE


HARD	Pre-preg		 <p>STM 22</p>	Carbon		 <p>STM 24 (1.5-2 mm)</p>	Simona		 <p>STM 25 (3 mm)</p>
	Polypropylene transparent		 <p>STM 23 (4 mm)</p>	Rhenoflex		 <p>STM 18 (1,8 mm)</p>	Metal toe cap		



HEEL COUNTER



TOE CAP





BRANDSOLE

EYELET · HOOK · TEXTILE HOOK · STUD · BUCKLE








QUICK LASING SYSTEM







WATER-REPELLENT ZIPPER H₂O ZIPP



EVA



RUBBER




LEATHER







EVA
09 White




EVA
19 Light beige



EVA
17 Dark beige



EVA
41 Taupe




EVA
56 Grey



EVA
46 Brown



EVA
81 Black



EVA
2 mm

RECENTLY ADDED

BR09-09	BR08-19	BR07-17	BR06-41	BR05-46	BR04-81
---------	---------	---------	---------	---------	---------

TPU
1,2 mm

BR23-09/82	BR24-09/78	BR33-09/81	BR25-56/81
------------	------------	------------	------------

TPU
DIAMOND
1,8 mm

BR29-09	BR32-81
---------	---------

TPU
LAVERO
3 mm

BR20-09	BR18-19	BR21-05
---------	---------	---------

BR19-17	BR15-56	BR17-41	BR16-46	BR14-81
---------	---------	---------	---------	---------

RUBBER
1,8 mm

BR31-87	BR30-81
---------	---------

RUBBER
STONE
Recycled
3 mm

BR27-81



SPORT



BR26-81




Super New Flex
35 Sh

- Durability 
- Lightness/softness 
- Anti-slip 



Blue Clasik 701
50 Sh

- Durability 
- Lightness/softness 
- Anti-slip 




Vibram 8870
45 Sh

- Durability 
- Lightness/softness 
- Anti-slip 




Astro-Star
50 Sh

- Durability 
- Lightness/softness 
- Anti-slip 



Swallow
45 Sh

- Durability 
- Lightness/softness 
- Anti-slip 



Vibram Mandorlo
45 Sh

- Durability 
- Lightness/softness 
- Anti-slip 



Astrolight
50 Sh

- Durability 
- Lightness/softness 
- Anti-slip 



Pluriball New Flex
50 Sh


- Durability 
- Lightness/softness 
- Anti-slip 



Stone Recycled
50 Sh



Contol
65 Sh



Topy Cellulux
55 Sh



Nylon Topy




Crepina Soft
60 Sh



Vibram Claw
60 Sh



Autoband
65 Sh



Resistant to oil and acid



Recycled material

Vibram 446



81

Vibram 1220



81

Vibram 1319



81

Vibram 2002



09

81

Vibram 2609



46

81

Vibram 2644



46

81

Vibram 2667



46

81

Vibram 2668



46

81

Vibram 2690



81

Vibram 528

- Durability
- Lightness/softness
- Anti-slip



46 81

Vibram 2603

- Durability
- Lightness/softness
- Anti-slip



81

Vibram 2655

- Durability
- Lightness/softness
- Anti-slip



46 81

Vibram 2678

- Durability
- Lightness/softness
- Anti-slip



81

Vibram 2810

- Durability
- Lightness/softness
- Anti-slip



46 81

Flensburg

- Durability
- Lightness/softness
- Anti-slip



81

Grip Line 505

- Durability
- Lightness/softness
- Anti-slip



46 81

Kiel

- Durability
- Lightness/softness
- Anti-slip



81

Lubeck

- Durability
- Lightness/softness
- Anti-slip



81



Resistant to oil and acid



Antistatic

Mammoth

- Durability: 5 orange squares
- Lightness/softness: 4 blue squares, 1 white square
- Anti-slip: 5 green squares




Sportflex

- Durability: 4 orange squares, 1 white square
- Lightness/softness: 4 blue squares, 1 white square
- Anti-slip: 5 green squares



5001

- Durability: 5 orange squares
- Lightness/softness: 4 blue squares, 1 white square
- Anti-slip: 5 green squares



Venus

- Durability: 4 orange squares, 1 white square
- Lightness/softness: 4 blue squares, 1 white square
- Anti-slip: 5 green squares




Ventosa

- Durability: 4 orange squares, 1 white square
- Lightness/softness: 4 blue squares, 1 white square
- Anti-slip: 5 green squares



Carrion

- Durability: 4 orange squares, 1 white square
- Lightness/softness: 4 blue squares, 1 white square
- Anti-slip: 5 green squares




Vibram 586

- Durability: 4 orange squares, 1 white square
- Lightness/softness: 5 blue squares
- Anti-slip: 5 green squares



Vibram 9105

- Durability: 5 orange squares
- Lightness/softness: 4 blue squares, 1 white square
- Anti-slip: 5 green squares



Roccia Block

- Durability: 5 orange squares
- Lightness/softness: 4 blue squares, 1 white square
- Anti-slip: 5 green squares



Resistant to oil and acid

Sole leather




SLM 02 (4-4,5 mm)


Anti-slip Sole - 52



Anti-slip Heel - 52



Anti-slip Sole




WOODEN AND LEATHER HEELS

Nr. 2



Nr. 6



LEFA



← Heel

← Heel tap

Shell PU1662

Durability

Lightness/softness

Anti-slip

Size
40-45



81

Nadia

Durability

Lightness/softness

Anti-slip

Size
36-41



81

Nel

Durability

Lightness/softness

Anti-slip

Size
37-41



81

Boston

Durability

Lightness/softness

Anti-slip

Size
39-45



81

Shore Hardness	Model	Dimensions	Material	Model	Dimensions	Material	Model	Dimensions	Material	
15-30 Sh	IM 24	4 mm 12 mm	Plastazote (15 Sh)	IM 23	6 mm	Multiform (15 Sh)	IM 21	1.6 mm	PPT (20 Sh)	
	IM 02	3.2 mm	PPT (20 Sh)	IM 26	6 mm	New flex (20 Sh)	IM 28	2 mm	Multiform (30 Sh)	
	IM 22	6.4 mm	PPT (20 Sh)	IM 29	1.5 mm	Poron (20 Sh)	IM 30	3 mm	Poron (20 Sh)	
	IM 13	15 mm 25 mm	Multiform (25 Sh)	IM 09	4 mm 6 mm	Nora (30 Sh)				
	35 Sh	IM 11	4 mm	Multiform	IM 14	5 mm 7.5 mm 10 mm 15 mm	Multiform	IM 27	5 mm 10 mm 20 mm	Multiform
		IM 06	5 mm 10 mm	Microcork	IM 31		EVA			

45 Sh

IM 32



55 Sh

IM 05

4 mm
6 mm
10 mm




IM 25

3 mm
5 mm
10 mm
15 mm



IM 33



BM04



For patients with diabetes and rheumatism

BM09



For patients with diabetes and rheumatism and cerebral paralysis

BM10



For patients with diabetes and rheumatism and cerebral paralysis

BM11



For daily orthopedic adult men insoles with weight over 75 kg

BM13



For daily orthopedic adult with extra midfoot support

NATURAL LEATHER

LN02

Calflining cream,
0.7-0.9 mm



LN45

Natura lining burro,
1.0-1.1 mm



LN26

Lamnappa naturel,
0.7-0.9 mm



LN22

Calflining camel,
0.9-1.1 mm



LN24

Calflining nero,
0.7-0.9 mm



LN32

Horse leather,
0.9-1.1 mm



MICROFABRIC

LN18

Suetex beige,
0.7 mm



LN17

Suetex grey,
0.7 mm



LN19

Micro suetex black,
0.6 mm



LN33

ON STEAM
perforated
beige,
0.8 mm



LN34

ON STEAM
perforated
black,
0.8 mm



THERMO

LN28

Natural fur
off white



LN30

Natural fur
brown



LN27

Natural fur
black



LN44

Foam
laminated
velour dark
grey,
4 mm



PADDED

LN50

Blue textile +
foam,
2.5 mm



LN51

ON STEAM
black + EVA
2.8 mm



LN53

ON STEAM
beige + poron
2.8 mm



BREATHABLE

LN47

DuraFit
PLUS,
1.1 mm



RECYCLED

LN46

Recycled
camel,
1 mm



Recycled material



ORTOPEDIJOS **TECHNIKA**

Krėvės pr. 129, LT-50312 Kaunas, Lithuania
E-mail: info@orthopedic-pro.com

Phone +370 37 312 686
Mobile +370 665 28 542

www.orthopedic-pro.com

