

FIXALTO®



ORTHOPEDIC ORTHOSES

2025

RehaFIX®

SpecialFIX™

SoftFIX®

BasicFIX™

Content

CERVICAL ORTHOSES

KR 0-12-1
BasicFIX™ Cervi Medium 3

SHOULDER ORTHOSES

RT 4-3-0
RehaFIX® Omo Light THERMO 3
 AP 2
RehaFIX® Omo Medium ABDUCTUS 3

THORACO-LUMBO-SACRAL ORTHOSES

KR 1-14-1
BasicFIX™ Thoraco Light RECLINE 8
 KR 1-14
BasicFIX™ Thoraco Medium RECLINE 8
 KR 1-17
BasicFIX™ Th-Lumbo Medium RECLINE 9
 KR 1-6
SpecialFIX™ Clavi Strong 9

ANKLE-FOOT ORTHOSES

SoftFIX® Malleo Light 2.0 17
 KT 0-2-2
BasicFIX™ Malleo Light AIR 17
 KT 0-4-1
BasicFIX™ Malleo Light THERMO 18
 KT 0-2-1
RehaFIX® Malleo Strong 18
 EASY – ON Slipper 18

SHOULDER-ELBOW ORTHOSES

RL 2-2
RehaFIX® Arm Slim Medium AIR 4
 RL 4
RehaFIX® Omo Strong DESAULT 4
 RT 3-5-4
SoftFIX® Cubi Light 2.0 4
 RT 3-5-1
BasicFIX™ Cubi Medium THERMO PLUS 5
SpecialFIX™ Epi Strap 5

LUMBO-SACRAL/LUMBAR ORTHOSES

KR 0-17-6
RehaFIX® Th-Lumbo Light 9
 KR 0-18-2
BasicFIX™ Lumbo Medium 10
 KR 0-17-3
SoftFIX® Lumbo Medium 10
 KR 0-20
SpecialFIX™ L-Sacro Strong 10
 KR 0-23
BasicFIX™ L-Sacro Strong 11
 KR 0-20-1
BasicFIX™ L-Sacro Strong THERMO 11
 KR 0-25
SpecialFIX™ Lumbo Strong PLUS 12
 KR 0-26
RehaFIX® L-Sacro Extra strong 12

MATERNITY ORTHOSES

KVP 6-6-1
SoftFIX® Lumbo Medium MAMA 19
 KVP 6-6
SoftFIX® Lumbo Strong MAMA 19
 FP 2
SpecialFIX™ Illium Strong ABDUCTUS 19

FINGER/WRIST/HAND ORTHOSES

RT 0-1-1
BasicFIX™ Digiti Light 5
 RT 0-1-2
BasicFIX™ Carpus Light 6
 RT 0-1-5
SpecialFIX™ Carpus Light THERMO 6
 RT 1-8-3
BasicFIX™ Carpus Light 2.0 6
 RT 1-8-6
SoftFIX® Manu Light 7
 RT 1-8-2
BasicFIX™ Manu Medium 7
 RT 1-8-4
RehaFIX® Manu Medium 7
 RT 1-8-8
RehaFIX® Manu Medium R/L 8

KNEE ORTHOSES

KT 3-8-14
BasicFIX™ Genu Light 13
SoftFIX® Genu Light 2.0 13
 KT 3-8-9
SpecialFIX™ Patella Medium THERMO SUBLUX 14
 KT 3-8-5
SpecialFIX™ Genu Medium THERMO 14
SpecialFIX™ Patella Strap 14
 KA 3-4 SER
SpecialFIX™ Genu Strong 15
 KA 3-3-1-2
RehaFIX® Genu Strong THERMO 15
 KA 3-3-1
BasicFIX™ Genu Strong THERMO 16
 KT 3-5-0
RehaFIX® Genu Extra strong AIR 16

OTHER ORTHOSES

SpecialFIX™ Stoma Belt 2.0 20
 BB 1-6
SpecialFIX™ Hernia Medium 20
 KR 0-22
SpecialFIX™ Sacro Strong 20

CUSTOM MADE ORTHOSES

Face mask 21
 Protective helmet 21

FIXALTO®

RehaFIX®

RehaFIX® stands for one of the groups of orthopaedic orthoses under the registered trademark FIXALTO of AB Ortopedijos technika – orthopaedic orthoses for the rehabilitation period after injuries or operations.

Reha means **rehabilitation, restoration, recovery** (movement, organ functions, etc.).
FIX denotes the function of the orthosis – to fix, stabilize the areas of the body under physical exertion and/or pain.

SpecialFIX™

SpecialFIX™ stands for one of the groups of orthopaedic orthoses under the registered trademark FIXALTO of AB Ortopedijos technika – orthopaedic orthoses with a special purpose and complex construction.

Special means **dedicated, extraordinary**.
FIX denotes the function of the orthosis – to fix, stabilize the areas of the body under physical exertion and/or pain.

SoftFIX®

SoftFIX® stands for one of the groups of orthopaedic orthoses under the registered trademark FIXALTO of AB Ortopedijos technika – orthopaedic 3D knitted orthoses for light fixation.

Soft means **gentle, light**.
FIX denotes the function of the orthosis – to fix, stabilize the areas of the body under physical exertion and/or pain.

BasicFIX™

BasicFIX™ stands for one of the groups of orthopaedic orthoses under the registered trademark FIXALTO of AB Ortopedijos technika – orthopaedic orthoses for prophylactic use with easy construction.

Basic means **basic, elementary**.
FIX denotes the function of the orthosis – to fix, stabilize the areas of the body under physical exertion and/or pain.

NECK



XXXXX – medium

Cervical orthosis KR 0-12-1

BasicFIX™ Cervi Medium

NEW DESIGN

Recommended for use:

- during post-surgery recovery period;
- during the treatment of neck area osteochondrosis, tendon, muscle and nerve diseases;
- during rehabilitation period after neck area spinal column traumas, strain and tear of muscles, ligaments and tendons, after vertebra movement and surgeries of neck area;
- after trauma (strain of neck ligaments).

SIZE	XS	S	M	L	XL	2XL
CM	30-32	32-34	34-37	37-41	41-44	44-47
HEIGHT	7,5	7,5	8,5	8,5	8,5	8,5



SHOULDER



XXXXX – light

Shoulder orthosis RT 4-3-0

RehaFIX® Omo Light THERMO

Recommended for use:

- during the treatment of tendon and ligament strains, bruising and inflammation of the shoulder joint;
- an orthosis stabilizes and heats up the shoulder joint, reduces strain and pain.

Product features:

made from air-tight fabric.

Note: the thickest part of humerus is measured.

SIZE	XS	S	M	L	XL	2XL	SIDE:
CM	24-26	26-29	29-32	32-35	35-38	38-42	UNIVERSAL



XXXXX – medium

Shoulder orthosis AP 2

RehaFIX® Omo Medium ABDUCTUS

Immobilizes arm in proper position.

Recommended for use:

- in case of humerus dislocation, shoulder joint bruising and inflammations;
- during post-surgery recovery period (after reconstructive surgeries of shoulder arch ligaments and muscles).

Product features:

- supports the shoulder joint in a physiologically correct position with external rotation of the shoulder;
- the product shape and stripes for position adjustment depending on the clinical situation;
- two bands for forearm fixation.

ONE SIZE | SIDE: UNIVERSAL

Shoulder - elbow orthosis RL 2-2

RehaFIX® Arm Slim Medium AIR



Immobilizes arm in proper position.

Recommended for use:

- after traumas and/or surgeries of shoulder bones, ligaments and muscles;
- in case of paralysis, inflammation of nerves, muscles, tendons and shoulder joint;
- during recovery period after shoulder joint dislocation (depending on the type of dislocation).

Product features:

the base is made of porous breathable fabric.



XXXXX – medium

	SIZE	XXS	XS	S	M	SIDE:
	CM	25-30	30-35	35-40	40-45	UNIVERSAL

Shoulder - elbow orthosis RL 4

RehaFIX® Omo Strong DESAULT



Immobilizes arm in proper position.

Recommended for use:

- after shoulder joint bruising or strain;
- after arthroscopic surgeries;
- in case of shoulder joint arthritis;
- after ruptures of collarbone and humerus proximal ends (without displacement).

Product feature:

the additional band for stronger fixation can be used.



XXXXX – strong

	SIZE	XXS	XS	S	M	SIDE:
	CM	< 30	30-35	35-40	> 40	UNIVERSAL

ELBOW

Elbow orthosis RT 3-5-4

SoftFIX® Cubi Light 2.0



The orthosis is anatomically shaped using modern 3D knitting technology. It fits snugly, easily and optimally presses and widely reinforces the elbow joint. The orthosis with unique honeycomb structure is highly breathable, elastic and functional.

Recommended for use:

in cases of epicondylitis, swelling caused by arthrosis/arthritis, after traumas and surgeries as well as for trauma prevention.

Product features:

- silicone inserts, attached in the inner side of the orthosis with seamless technology, help to evenly distribute the load on the surrounding soft tissues;
- the orthosis will always remain in the proper place (will not slip off) as these inserts repeat the anatomy of joint and fix the orthosis on humerus condyles.



XXXXX – light

Note: measure the thickest part of the forearm.

	SIZE	S	M	L	XL	2XL	SIDE:	COLOUR:
	CM	21-23	23-25	25-27	27-29	29-31	UNIVERSAL	GREY OR BEIGE

ELBOW



XXXXX – medium

Elbow orthosis RT 3-5-1

BasicFIX™ Cubi Medium THERMO PLUS

Compensates the unrestraint of elbow ligaments. Provides warmth to the elbow joint.

Recommended for use:

- instability of elbow joint;
- inflammation of muscles and tendons;
- bursitis.

Product features:

with immobilizing straps.

Note: measure the thickest part of the forearm.



SIZE	XS	S	M	L	XL	SIDE:
CM	23-25	25-27	27-30	30-33	33-36	UNIVERSAL



XXXXX – strong

Elbow orthosis

SpecialFIX™ Epi Strap

Partly immobilizes the area of tendons attaching to the external humeral condyle.

Recommended in case of external or internal epicondylitis.

Especially useful for tennis players and masseurs.

For forearm circumferences from 18-36 cm.

Note: measure the thickest part of the forearm.



ONE SIZE | SIDE: UNIVERSAL



FINGER



XXXXX – light

Finger orthosis RT 0-1-1

BasicFIX™ Digiti Light

Recommended for use:

- after first metacarpal bone traumas (ruptures, bruises, strains of finger ligaments);
- in case of tendon inflammation.

Product feature:

plastic fastener for motion limitation.



SIZE	XS	S	M	L	SIDE:
CM	12-14	14-16	16-18	18-20	LEFT OR RIGHT

Wrist orthosis RT 0-1-2

BasicFIX™ Carpus Light

Elastic wrist orthosis with contact closure straps optimally stabilizes the wrist in the right position.

Recommended for use:

- to protect the wrist joint from the big load and physical exertion;
- during rehabilitation period after ligament and tendon strains and tear;
- in case of tendinitis, wrist joint arthrosis.

For wrist circumferences from 14-28 cm.



ONE SIZE

SIDE: LEFT OR RIGHT



XXXXX - light

Wrist orthosis RT 0-1-5

SpecialFIX™ Carpus Light THERMO

Easily immobilizes and warms wrist joint.

Recommended for use:

- in case of wrist instability, wrist joint arthrosis, after wrist joint ligament strain,
- during rehabilitation period after wrist joint canal surgery.

For wrist circumferences from 14-28 cm.



ONE SIZE

SIDE: LEFT OR RIGHT



XXXXX - light

Wrist orthosis RT 1-8-3

BasicFIX™ Carpus Light 2.0

Immobilizes and warms wrist joint.

Recommended for use:

- in case of wrist joint canal syndrome;
- after strain of ligaments and tendons and during other cases.

For wrist circumferences from 13-25 cm.



ONE SIZE

SIDE: UNIVERSAL



XXXXX - light



XXXXX – light

Wrist - hand orthosis RT 1-8-6

SoftFIX® Manu Light

The orthosis is anatomically shaped with silicon insert using modern 3D knitting technology. It fits snugly, easily and optimally presses and lightly reinforces the wrist and the hand. The orthosis with unique honeycomb structure is highly breathable, elastic and functional.

Recommended for use:

- during rehabilitation period after wrist and lower forearm traumas (ruptures and ligament strains);
- in case of wrist carpal tunnel syndrome;
- in case of wrist tendinitis, tendovaginitis;
- after lower forearm or wrist fracture.



Note. For easier fitting of this product, we recommend to order the Easy-ON Slipper (see page 18).

Product features:

- silicone insert, attached in the inner side of the orthosis with seamless technology, helps to evenly distribute the load on the surrounding soft tissues;
- additional elastic strap for fixation strength control.



SIZE	S	M	L	SIDE:
CM	14-16	16-18	18-20	LEFT OR RIGHT



XXXX – medium

Wrist - hand orthosis RT 1-8-2

BasicFIX™ Manu Medium

Immobilizes wrist and hand in proper position.

Recommended for use:

- after wrist and forearm traumas (ruptures and ligament strains);
- in case of forearm tendon inflammation.

Product features:

- rigid metal splint and 2 plastic side fasteners for limitation of motion;
- double adjustable fastener helps to evenly distribute the load on the wrist and adjust the fixation strength.



SIZE	XS	S	M	L	SIDE:
CM	12-14	14-16	16-18	18-20	LEFT OR RIGHT



XXXX – medium

Wrist - hand orthosis RT 1-8-4

RehaFIX® Manu Medium

Supports the wrist and palm in a physiologically correct position after injuries and in wrist and forearm pathologies.

Recommended for use:

- in case of tendovaginitis, strain and tear of wrist ligaments;
- after knucklebone ruptures;
- after spoke-bone rupture in typical location and wrist area.

Product features:

- rigid metal splint and 2 plastic side bones for limitation of motion;
- triple adjustable fastener helps to evenly distribute the load on the wrist and adjust the fixation strength;



SIZE	XS	S	M	L	SIDE:
CM	12-14	14-16	16-18	18-20	LEFT OR RIGHT

Wrist - hand orthosis RT 1-8-8
RehaFIX® Manu Medium L/R



Supports the wrist and palm in a physiologically correct position after injuries and in wrist and forearm pathologies.

Recommended for use:

- in case of tendovaginitis, strain and tear of wrist ligaments;
- after knucklebone ruptures;
- after spoke-bone rupture in typical location and wrist area.

Product features:

- rigid metal splint and 2 plastic side bones for limitation of motion;
- triple adjustable fastener helps to evenly distribute the load on the wrist and adjust the fixation strength;



XXXXX - medium



SIZE	XS	S	M	L	SIDE:
CM	12-14	14-16	16-18	18-20	UNIVERSAL

Thoraco - lumbo - sacral orthosis KR 1-14-1

BasicFIX™ Thoraco Light RECLINE



Recommended for use:

- to correct defective postures;
- in conjunction with other therapies;
- in the treatment of spinal columns in the thorax area.



SIZE	XXS	XS	S	M	L
CM	55-65	65-75	75-85	85-95	95-105



XXXXX - light

Thoraco - lumbo - sacral orthosis KR 1-14

BasicFIX™ Thoraco Medium RECLINE




Recommended for faulty posture correction.

Product features:

- stainless steel bones evenly distribute the load on the spine;
- elastic straps help to stretch and maintain proper posture.



SIZE	XS	S	M	L	XL	2XL	3XL	4XL
a, CM	50-60	60-70	70-80	80-90	90-100	100-110	110-120	120-130
h, CM	36-41	38-43	40-45	42-47	44-49	46-51	47-52	48-53



XXXXX - medium



XXXXX – medium

Thoraco - lumbo - sacral orthosis KR 1-17
BasicFIX™ Th-Lumbo *Medium* RECLINE





Recommended for use:

- to correct faulty posture;
- in case of juvenile osteochondrosis of spine;
- in case of minor kyphosis, lordosis and scoliosis;
- during recovery period after light thoracic spine area traumas, bruising and strain;
- in case of inflammation of nerves, muscles and ligaments.

Product features:

- stainless steel bones evenly distribute the load on the spine;
- elastic straps help to stretch and maintain proper posture.

SIZE	XS	S	M	L	XL	2XL	3XL	4XL
a, CM	50-60	60-70	70-80	80-90	90-100	100-110	110-120	120-130
h, CM	36-41	38-43	40-45	42-47	44-49	46-51	47-52	48-53



XXXXX – strong

Cervical - breast orthosis KR 1-6
SpecialFIX™ Clavi *Strong*



Recommended for use:

- to stabilize clavicle – improves posture;
- in case of small deformations of the spinal column;
- in case of various deformations in the breast area;
- after fracture of the collarbone.

SIZE	XXS	XS	S	M	L	XL	2XL
CM	25-30	30-35	35-40	40-45	45-50	50-55	55-60



XXXXX – light

Lumbo - sacral orthosis KR 0-17-6
RehaFIX® Th-Lumbo *Light*

Recommended for use:

- after stomach operations, or for muscle splitting and muscle weakness after pregnancy;
- to prevent ruptures in the abdominal area.

SIZE	S	M	L
CM	60-80	80-100	100-120



Lumbo - sacral orthosis KR 0-18-2 (with double fastening)

BasicFIX™ Lumbo Medium


Recommended for use:

- in case of pain in lumbar section and radiculitis;
- in case of strain of hip section muscles and ligaments;
- after abdominal surgeries, muscle rupture, in case of abdominal weakness after pregnancy and to support abdominal hernia.

Product features:

- double fastener for fixation strength adjustment;
- 2 metal fasteners to stabilize the lumbar spine area and evenly distribute the load on the spine.

Note: Choose one size smaller orthosis immediately after birth than recommended in the size chart.



XXXXX – medium



SIZE	XS	S	M	L	XL	2XL	3XL	4XL
CM	60-70	70-80	80-90	90-100	100-110	110-120	120-130	130-140

Lumbo - sacral orthosis KR 0-17-3

SoftFIX® Lumbo Medium

The orthosis is anatomically shaped with silicon pad using modern 3D knitting technology. It fits snugly, optimally presses and reinforces the lumbar part of the spine. The orthosis with unique honeycomb structure is highly breathable, elastic and functional.

Product features:

- 2 plastic fasteners fitted on the inside evenly distribute the load on the spine;
- special removable silicone spiky pad massages the lumbar – sacral area and stimulates recovery processes;
- seamless technology guarantees wearing comfort.

Recommended for use:

in cases of pain in the lumbar section, radiculitis, strain of muscles and ligaments in the lumbar section as well as to protect the lumbar section.



XXXXX – medium



SIZE	S	M	L	XL	2XL
CM	70-80	80-90	90-100	100-110	110-120



Lumbo - sacral orthosis KR 0-20

SpecialFIX™ L-Sacro Strong

Recommended for use:

- after lumbar section traumas (bruises, muscle strains);
- in case of diseases of intervertebral disks of lumbar spine area (spondilolistesis, osteochondrosis);
- during rehabilitation period after non-complicated rupture of vertebral column;
- after abdominal and renal surgeries.

Product features:

- additional stripes for fixation strength adjustment;
- 2 stainless steel bones on the back and 2 spiral steel bones on the sides to stabilize the lumbar spine area and evenly distribute the load on the spine.



XXXXX – strong



SIZE	XS	S	M	L	XL	2XL	3XL
CM	60-70	70-80	80-90	90-100	100-110	110-120	120-130



XXXX – strong

Lumbo - sacral orthosis KR 0-23

BasicFIX™ L-Sacro Strong

Recommended for use:

- in cases of pain in the lumbar section, radiculitis, strain of muscles and ligaments in the lumbar section (spondilolistesis, osteochondrosis);
- after abdominal surgeries, muscle rupture, in case of abdominal weakness after pregnancy and to support abdominal hernia;
- in case of intervertebral disk diseases;
- to treat spina bifida.

Product features:

- double fastener for fixation strength adjustment;
- cross-bands for additional fixation and lumbar spine support;
- 2 stainless steel and 2 plastic bones in the back-to stabilize the lumbar spine area and evenly distribute the load on the spine.



SIZE	XS	S	M	L	XL	2XL	3XL
CM	60-70	70-80	80-90	90-100	100-110	110-120	120-130



XXXX – strong

Lumbo - sacral orthosis KR 0-20-1

BasicFIX™ L-Sacro Strong THERMO

Recommended for use:

- in case of pain in lumbar section and radiculitis;
- in case of strain of lumbar section muscles and ligaments;
- to protect lumbar section.

Product features:

- made of material, featuring compression and warming effect;
- additional detachable stripes for enhanced lumbar stability and fixation strength adjustment;
- 4 plastic bones in the back- to stabilize the lumbar spine area and evenly distribute the load on the spine.



SIZE	XS	S	M	L	XL	2XL
CM	60-70	70-80	80-90	90-100	100-110	110-120

Lumbo - sacral orthosis KR 0-25 (with extra shoulder-strap)

SpecialFIX™ Lumbo Strong PLUS

The orthosis provides average fixation of the lumbar spine. It stabilizes ligament-articular apparatus. The product eliminates imbalance of paravertebral muscles in the lumbar spine and reduces the load on the lumbar muscles.

Recommended for use:

- in case of lumbar pain due to muscle imbalance, scoliosis, shortening of one of the lower extremities, radiculitis, osteochondrosis, spondylarthrosis, protrusions and hernia of intervertebral disks, dislocation of vertebra;
- during rehabilitation after lumbar spine injuries and operations;
- in case of splanchnoptosis of abdominal cavity;
- during functional therapy in combination with manual therapy, massage, medication;
- for the prevention of injuries and diseases of waist in case of long-time loads;
- to maintain therapeutic effect achieved during the course of manual therapy and therapeutic exercises.



XXXXX – strong

Product features:

the product is equipped with metal modulated stiffening plates and additional adjustable tightening.



SIZE	XS	S	M	L	XL	2XL
(a+b)/2, CM	<70	70-85	85-97	97-105	105-115	>115



Lumbo - sacral orthosis KR 0-26

RehaFIX® L-Sacro Extra strong

This product provides strong fixation of the lumbosacral spine. It restricts sudden movements.

Recommended for use:

- after injuries and operations;
- in case of lumbar pain due to radiculitis, sharpening of osteochondrosis;
- to maintain therapeutic effect achieved during the course of manual therapy and therapeutic exercise;
- to prevent complications in treatment of lumbar spine diseases.



XXXXX – extra strong

Product features:

the product is equipped with metal modulated stiffening plates and additional adjustable tightening.



SIZE	XS	S	M	L	XL	2XL	HEIGHT
(a+b)/2, CM	<70	70-85	85-97	97-105	105-115	>115	>165





XXXXX - light

Knee orthosis KT 3-8-14

BasicFIX™ Genu Light

The orthosis stabilizes knee-cap and warms knee joint.

Recommended for use:

in cases of strain, inflammation of knee-cap and knee joint ligaments, tendons and muscles.

For knee circumferences from 28-48 cm.



ONE SIZE

SIDE: UNIVERSAL



XXXXX - light

Knee orthosis

SoftFIX® Genu Light 2.0



The orthosis is anatomically shaped using modern 3D knitting technology. It fits snugly, easily and optimally presses and widely reinforces the knee joint. The orthosis with unique honeycomb structure is highly breathable, elastic and functional.

Recommended for use:

in cases of knee pain, diseases of joint cartilages, faulty position of the knee patella, strain of the knee joint ligaments and muscles, arthrosis, osteoarthritis, sinovitis, chondritis, patella osteochondropathies of various origins, inflammation of tendons and joint bursae of knee joint section, osteochondropathy of tibial tuberosity, instability of patella, inflammation in the patellar ligament as well as for trauma prevention.

Product features:

- spiral steel bones secured on the inside with seamless technology for enhanced support;
- silicone ring stabilizes the kneecap in a correct position; because of the silicone ring the brace will always remain in place. Knee joint load is uniformly distributed.

Note: circumference measured 14 cm above and 12 cm below the middle of the kneecap with a given flexion of 30°.

SIZE	S	M	L	XL	2XL	SIDE:	COLOUR:
a, CM (above)	41-44	44-47	47-50	50-53	53-56	UNIVERSAL	GREY OR BEIGE
b, CM (below)	31-34	34-37	37-40	40-43	43-46		



Knee orthosis KT 3-8-9

SpecialFIX™ Patella Medium THERMO SUBLUX

Recommended for use:

- in case of unrestraint, strain, tear, partial or full dislocation of patella ligaments;
- during recovery period after plastic surgeries of patella ligaments.

Product features:

lateral support on the external side of the knee stabilizes the patella in a correct position.

Support inside is covered with natural leather.



XXXXX – medium



SIZE	XS	S	M	L	XL	SIDE:
CM	30-33	33-36	36-40	40-44	44-46	UNIVERSAL

Knee orthosis KT 3-8-5

SpecialFIX™ Genu Medium THERMO

Stabilizes knee joint and patella.

Recommended for use:

in case of strain and tear knee joint ligaments.

Product features:

- made of material, featuring compression and warming effect;
- 4 spiral steel bones (2 on each side) for enhanced stability;
- additional detachable straps for fixation strength adjustment.



XXXXX – medium



SIZE	XS	S	M	L	XL	2XL	SIDE:
CM	30-33	33-36	36-40	40-44	44-46	46-50	UNIVERSAL

Knee orthosis

SpecialFIX™ Patella Strap

Equally allocates the load directed to patella ligaments.

Recommended in case of:

- osteochondropathy of tibial tuberosity (Osgood-Schlatter disease);
- strain of patella ligaments.

Product feature:

support inside is covered with natural leather.

For kneecap circumferences from 25-43 cm.



XXXXX – strong



ONE SIZE	SIDE: UNIVERSAL
----------	-----------------



XXXX – strong

Knee orthosis KA 3-4 SER SpecialFIX™ Genu Strong



Stabilizes knee joint at the required angle or to limit the range of motion.

Recommended for use:

- in case after ligament tears, in case of spastic cerebral palsy, hemiplegia, paraplegia and tetraplegia;
- during rehabilitation period after traumas and operations.

Product features:

- rigid metal hinges on the sides for secure stability;
- the orthosis has an angle adjustment function;
- open type: easy to put on and take off as the orthosis opens out fully.



SIZE	S	M	L	XL	2XL	3XL	SIDE:
CM	33-36	36-40	40-44	44-46	46-50	50-55	UNIVERSAL



XXXX – strong

Knee orthosis KA 3-3-1-2 RehaFIX® Genu Strong THERMO

Immobilizes the knee joint and centres the patella. Firm lateral splints particularly fasten the knee joint.

Recommended in case of:

- damage (strain) of collateral ligaments;
- gonarthrosis, osteochondropathy (Osgood-Schlatter disease);
- various forms of cerebral palsy (tetraplegia, paraplegia).

Product features:

- made of material featuring compression and warming effect;
- rigid metal hinges on the sides for secure stability;
- silicone ring stabilizes the kneecap in a correct position; because of the silicone ring the orthosis will always remain in place. Knee joint load is distributed uniformly;
- additional detachable stripes for fixation strength adjustment;
- open type: easy to put on and take off as the orthosis opens out fully.



SIZE	XS	S	M	L	XL	2XL	SIDE:
CM	30-33	33-36	36-40	40-44	44-46	46-50	UNIVERSAL

Knee orthosis KA 3-3-1

BasicFIX™ Genu Strong THERMO

Immobilizes the knee joint and centres the patella. Firm lateral splints particularly fasten the knee joint.

Recommended in case of:

- damage (strain) of collateral ligaments;
- gonarthrosis, osteochondropathy (Osgood-Schlatter disease);
- various forms of cerebral palsy (tetraplegia, paraplegia).

Product features:

- made of material featuring compression and warming effect;
- rigid metal hinges on the sides for secure stability;
- silicone ring stabilizes the kneecap in a correct position; because of the silicone ring the orthosis will always remain in place. Knee joint load is distributed uniformly;
- additional detachable stripes for fixation strength adjustment.



XXXXX – strong



SIZE	XS	S	M	L	XL	2XL	SIDE:
CM	30-33	33-36	36-40	40-44	44-46	46-50	UNIVERSAL

Knee orthosis KT 3-5-0

RehaFIX® Genu Extra strong AIR

Recommended for use:

- to provide first aid after knee joint traumas;
- during the treatment of inflammation, strain of ligaments and tendons, hemarthrosis and gonarthrosis;
- during recovery period after knee joint surgeries.

Product features:

- fastening straps and 4 inflexible metal splints (2 on the sides and 2 on the back) located inside stabilize knee joint. Therefore the orthosis can be worn instead of plaster splints;
- the base is made of porous breathable fabric.



XXXX – extra strong



SIZE	XS	S	M	L	XL	SIDE:
a, CM	38-43	44-48	49-53	54-58	59-63	UNIVERSAL
b, CM	28-32	33-36	37-40	41-44	45-48	



Ankle - foot orthosis SoftFIX® Malleo Light 2.0

The orthosis is anatomically shaped using modern 3D knitting technology. It fits snugly, easily and optimally presses and widely reinforces the ankle joint. The orthosis with unique honeycomb structure is highly breathable, elastic and functional.

Recommended for use:

in cases of swelling caused by arthrosis and arthritis, tendinomiopathy, week tarsus ligaments, damage and inflammation of the ankle joint ligaments, tendons and muscles as well as for trauma prevention.

Product feature:

silicone inserts fitted on the inside using seamless technology stabilize the ankle ligaments.



Note. The product comes with a slip for easier putting on.



SIZE	S	M	L	XL	SIDE:	COLOUR:
CM	19-21	21-23	23-25	25-27	UNIVERSAL	GREY OR BEIGE

XXXXX - light



Ankle - foot orthosis KT 0-2-2 BasicFIX™ Malleo Light AIR

Immobilizes the ankle joint.

Recommended for use:

- during recovery period after ankle joint traumas and surgeries;
- in case of inflammation of ankle joint or joint instability.

Product features:

- 4 plastic side bones (2 on each side)- for enhanced stability.
- the base is made of porous breathable fabric.



Note. For easier fitting of this product, we recommend to order the Easy-ON Slipper (see page 18).



SIZE	XS	S	M	L	XL	SIDE:
CM	20-24	24-28	28-32	32-36	36-39	UNIVERSAL

XXXXX - light

Ankle - foot orthosis KT 0-4-1

BasicFIX™ Malleo Light THERMO

Provides warmth to the ankle joint and Achilles tendon.

Recommended for use:

in case of damage and inflammation of ankle joint ligaments, tendons and muscles.

Product feature:

made of material featuring compression and warming effect;



Note. For easier fitting of this product, we recommend to order the Easy-ON Slipper (see page 18).



XXXXX – light



SIZE	XS	S	M	L	XL	SIDE:
CM	20-24	24-28	28-32	32-36	36-39	UNIVERSAL

Ankle - foot orthosis KT 0-2-1

RehaFIX® Malleo Strong

Recommended for use:

- in case of damages to ankle joint ligaments;
- tendons and muscles (bruising, strain, tear, joint instability).

Product features:

- 4 plastic side bones (2 on each side)- for enhanced stability;
- additional detachable straps for fixation strength adjustment.



SIZE	XS	S	M	L	XL	SIDE:
CM	20-24	24-28	28-32	32-36	36-39	UNIVERSAL



XXXXX – strong

EASY - ON slipper

(application for ankle - foot and wrist - hand supports)

Please put on the ankle or wrist support early in the morning after getting up if possible (or as necessary).

1. Slide your foot or hand into the provided easy – on slipper.
2. Then slide the ankle or wrist support over the donning aid, smoothing it up evenly.
3. Please make sure that the ankle or wrist support is not twisted and the seam is centered over the heel or hand.
4. The easy – on slipper can now be removed easily.
5. The product can be positioned as needed now.



ONE SIZE



XXXXX – medium

Lumbar orthosis KVP 6-6-1 (maternity belt)

SoftFIX® Lumbo *Medium* MAMA

Anatomically formed knitted orthosis gently provides support for the abdomen and reduces the pressure directed to the waist (back) during pregnancy; equally allocates the weight along the back and helps to avoid the back (spinal) and waist pains. Helps to sustain the correct set-up and to move actively. Protects from the muscles overstrain. Orthosis is elastic, breathable and comfortable to wear.

Recommended to be worn from the mid-pregnancy.



SIZE	S	M	L	XL	2XL
CM	90-100	100-110	110-120	120-130	130-140



XXXXX – strong

Lumbar orthosis KVP 6-6 (maternity belt)

SoftFIX® Lumbo *Strong* MAMA

Pregnancy belt is anatomically shaped using modern 3D knitting technology. It fits snugly and optimally holds the tummy during pregnancy. Helps to sustain correct posture and to facilitate the movement. Protects against muscle strain and formation of stretch marks. The orthosis with unique honeycomb structure is highly breathable, elastic and functional.

Recommended to be worn from the mid-pregnancy.

Product features:

- double hooking offers individual adjustment option;
- special rigidity reinforcements in the back area help to evenly distribute the waist (back) load and avoid lower back pain;
- seamless technology guarantees wearing comfort.



SIZE	S	M	L	XL	2XL
CM	90-100	100-110	110-120	120-130	130-140



XXXXX – strong

Hip abduction orthosis FP 2

SpecialFIX™ Illium *Strong* ABDUCTUS

Immobilizes hip joints in proper position.

Recommended for use:

in case of dysplasia in hip joint.

Note: product colours and material may be different from those shown in the picture.



SIZE	XXS	XS	S	M	L	XL	2XL	3XL	4XL	5XL
CM	14	16	18	20	22	24	26	28	30	32

Ostomy belt

SpecialFIX™ Stoma Belt 2.0



A belt with a unique design for patients with a stoma.

Recommended for use:

- the belt is made of elastic cotton fabric that does not irritate the skin and is reinforced in the stoma area with a mesh that limits elasticity where it is needed: the reinforced area presses the collector plate and prevents it from sticking to the skin;
- special pocket to hold the ostomy bag. The size of the pocket can be adjusted depending on the need, i.e. bag filling;
- the opening in the belt is adapted to various sizes of the openings of the stoma content collectors;
- the belt is suitable for a stoma originating on the left or right side, above or below the umbilicus.



SIZE	S/M	L/XL	SIDE:
CM	70-100	100-130	UNIVERSAL



Groin hernia truss BB 1-6

SpecialFIX™ Hernia Medium

Recommended for use:

- in case of the treatment of groin ruptures;
- to prevent ruptures from reoccurring;
- in both cases for ruptures that have not been operated.



SIZE	XS	S	M	L	SIDES:
CM	90-95	95-100	100-105	105-110	LEFT, RIGHT OR TWO - SIDED

XXXXX – medium

Sacro - iliac orthosis KR 0-22

SpecialFIX™ Sacro Strong

Immobilizes sacral spine and pelvic area.

Recommended for use:

- in case of rupture in symphysis pubis;
- during recovery period after surgeries in this area.

Product features:

- side fasteners for fixation strength adjustment;
- additional fasteners on the bottom always keep the orthosis in the right place.



SIZE	XS	S	M	L	XL	2XL
CM	85-90	90-95	95-100	100-105	105-110	110-115

XXXXX – strong



Face mask

Recommended for:

sportsmen after various facial injuries (nose or other facial bone fractures). With this protection, the athlete can return to the sport faster, which is extremely important for professionals.

Features of the protective face mask:

- made individually,
- protects against repeated injuries,
- extremely strong and light,
- the inner part of the mask is covered with a soft material that absorbs part of the impact and increases wearing comfort.

With this protection:

a professional athlete can resume playing in a matter of weeks instead of missing the rest of the season.

Production time:

3-5 working days.
If you need it earlier, please contact your export manager.

How to order?

Send us a 3D scan.

Protective carbon fiber face masks are not just for professional athletes, they can also help many active people who have experienced facial injury or surgery.



Protective helmet



Recommended for use:

patients who injure themselves or are unbalanced and may accidentally faint.

Product features:

- natural leather finishing;
- different helmet shapes (open, open front, open back or closed);
- three secure fastening options (plastic, metal or magnetic buckle);
- choice of the preferred colour from the palette in the order sheet.

Order flow without trial stage:

- helmets for fitting (3 sizes) – these are held by our client's company and are used only for measuring;
- protective helmet for wearing – individually made for the patient.

The advantage – no need to send the helmet for fitting, saving both time and transport costs.

How does it work?

1. Set of fitting helmets is sent to the client company.
2. To start the production of the custom-made helmet, we need to receive one of the three following options:
 - a photo of the patient wearing the fitting helmet;
 - measurements according to the order sheet;
 - the old helmet as a sample.
3. Patient receives a custom-made protective helmet to wear.


Important to know:

- the hair will be pinned when measuring.
- if fitting helmets are used, only photos of the patient wearing the fitting helmet are needed.
- if fitting helmets are not available, measurements or the old helmet can be sent instead. All options can also be provided together for the best result.



ORTOPEDIJOS **TECHNIKA**

*Export partners in 25 countries.
Cooperation with trust!*



Export department:

Krėvės pr. 129, LT-50312 Kaunas, Lithuania

Headquarters address:

Partizanų str. 17, LT-49476 Kaunas, Lithuania

E-mail: info@orthopedic-pro.com

Phone +370 37 312 686

Mobile +370 665 28 542

www.orthopedic-pro.com